



MICROLUX™


TRANSILLUMINATOR



Anterior Caries



Posterior Caries



Subgingival Calculus

The Microlux™ transilluminator provides a focused beam of cool white light for use in clinical diagnosis of caries, calculus, crown fractures and a host of other applications. The light is powered by a high intensity light emitting diode (LED) and transmitted through a focused glass fiber

optic element. By using the technique of fiber optic transillumination, the scope of oral diagnosis can be greatly improved without the introduction of additional x-rays. It also functions as an auxiliary light source which aids in general operative procedures, preventive dentistry and patient education.



AdDent, Inc. ■ 43 Miry Brook Road ■ Danbury, CT 06810-7414
Phone: (203) 778-0200 ■ Fax: (203) 792-2275 ■ www.addent.com

Anterior Caries



1a



1b




1c


Posterior Caries



2a



2b



2c

Subgingival Calculus



3a



3b

Transillumination Technique:

To detect anterior caries, place the probe on the labio-cervical region of the tooth and examine from the lingual with a mouth mirror, as shown in figure 1b. To diagnose posterior proximal caries, place the probe on the bucco-cervical area of the tooth. Light passes into the cervical region of the tooth structure and then radiates occlusally. Caries appear as dark triangular shadow on the occlusal surface as shown in figure 2b. Subgingival calculus can be demonstrated, as shown in figure 3b. It appears dark in contrast to the translucent tooth structure.

Features:

- White light emitting diode (LED) light source.
- Focused glass fiber optic light guide.
- Battery Operated.
- Measures only 5" x 5/8" diameter (13cm x 1.5cm diameter).

Benefits:

- Provides high intensity light for diagnostic procedures with extremely long 50,000-hour life and low battery drain.
- Autoclaveable light guide provides concentrated 3mm spot.
- Easily portable.
- Compact unit easily fits in your pocket.

References:

- Friedman, J. and Marcus, M.I.: Transillumination of the Oral Cavity with use of Fiber Optics. *J.A.D.A* 80: 801-809, April 1970.
- Cortes, D.F., Ekstrand, K.R. et al.: An In Vitro Comparison of the Ability of Fiber-Optic Transillumination, Visual Inspection and Radiographs to Detect Occlusal Caries and Evaluate Lesion Depth: *Caries Research* 2000, 34: 443-447.
- Davies, G.M., Worthington, H.V., et al: *British Dental Journal* 191: 145-147, August 2001.



MicroLux™ Transilluminator (Complete Kit) P/N 110019

- (1) - LED Transilluminator P/N 640006
- (1) - Glass Light Guide P/N 630003
- (3) - Batteries P/N 630006
- (1) - Pocket Clip P/N 640007



AdDent, Inc. ▪ 43 Miry Brook Road ▪ Danbury, CT 06810-7414
Phone: (203) 778-0200 ▪ Fax: (203) 792-2275 ▪ e-mail: mail@addent.com