









**WHAT
FOR
WHAT**

AIR-FLOW® POWDERS

EMS⁺

AIR - FLOW® POWDER		
ORIGINAL AIR - FLOW® THERAPY → SUPRAGINGIVAL	CLINICAL APPLICATIONS	AIR - FLOW® HANDPIECES
 <p>AIR - FLOW® POWDER CLASSIC</p> <ul style="list-style-type: none"> > Sodium bicarbonate > Average grain size ~65µm > Available in 6 different tastes <p>300 g bottle</p>	<ul style="list-style-type: none"> → REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS → SURFACE PREPARATION BEFORE DENTAL TREATMENTS 	 <p>AIR - FLOW®</p>
 <p>AIR - FLOW® POWDER CLASSIC COMFORT</p> <ul style="list-style-type: none"> > Sodium bicarbonate > Average grain size ~40µm <p>250 g bottle</p>	<ul style="list-style-type: none"> → REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS WITH OPTIMIZED PATIENT COMFORT → SURFACE PREPARATION BEFORE DENTAL TREATMENTS 	
 <p>AIR - FLOW® POWDER SOFT</p> <ul style="list-style-type: none"> > Glycine-based > Average grain size ~65µm <p>200 g bottle</p>	<ul style="list-style-type: none"> → REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS → TOOTH POLISHING 	
→ SUB - AND SUPRAGINGIVAL		
 <p>AIR - FLOW® POWDER PLUS</p> <ul style="list-style-type: none"> > Erythritol-based > Average grain size ~14µm > Contains chlorhexidine (0,3%) <p>100 g bottle</p>	<ul style="list-style-type: none"> → REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS WITH OPTIMIZED PATIENT COMFORT → TOOTH POLISHING <ul style="list-style-type: none"> → REMOVAL OF BIOFILM FROM POCKETS DOWN TO 4mm ON NATURAL TEETH AND IMPLANTS → REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, FROM 4 TO 10mm ON NATURAL TEETH AND IMPLANTS 	 <p>AIR - FLOW® AIR - FLOW® PLUS</p>
→ SUBGINGIVAL		
 <p>AIR - FLOW® POWDER PERIO</p> <ul style="list-style-type: none"> > Glycine-based > Average grain size ~25µm <p>120 g bottle</p>	<ul style="list-style-type: none"> → REMOVAL OF BIOFILM FROM POCKETS DOWN TO 4mm ON NATURAL TEETH AND IMPLANTS → REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, FROM 4 TO 10mm ON NATURAL TEETH AND IMPLANTS 	 <p>PERIO - FLOW® AIR - FLOW® PERIO</p>



NEW



NEW

**"I FEEL
GOOD"**

EMS-SWISSQUALITY.COM

EMS
ELECTRO MEDICAL SYSTEMS SA
Chemin de la Vuarpilliére 31
CH-1260 Nyon, Switzerland

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-dent.com